USER MANUAL LIST

KULLANIM KILAVUZU LİSTESİ

PRODUCT/ÜRÜN	PAGE/SAYFA
<i>NANO HİBRİT ÜNİVERSAL KOMPOZİT (tr)</i>	2
NANO HYBRID UNIVERSAL COMPOSITE(eng)	3
<i>KIDS NANO HİBRİT ÜNİVERSAL KOMPOZİT (tr)</i>	4
KIDS NANO HYBRID UNIVERSAL COMPOSITE (eng)	5
<i>NANO HİBRİT POSTERİOR KOMPOZİT (tr)</i>	6
NANO HYBRID POSTERIOR COMPOSITE (eng)	7
<i>AKIŞKAN (FLOW) KOMPOZİT (tr)</i>	8
FLOW COMPOSITE (eng)	9
<i>KIDS AKIŞKAN (FLOW) KOMPOZİT (tr)</i>	10
KIDS FLOW COMPOSITE (eng)	11
<i>BULK FLOW KOMPOZİT (tr)</i>	12
BULK FLOW COMPOSITE (eng)	13
<i>FLOR LINER KOMPOZİT (tr)</i>	14
FLOR LINER COMPOSITE (eng)	15
<i>ORTODONTİK AKIŞKAN (FLOW) ADEZİV (tr)</i>	16
ORTHODONTIC FLOW ADHESIVE (eng)	17
<i>BOND (7.JENERASYON DENTAL ADEZİV (tr)</i>	18
BOND (7th. GENER.DENTAL ADHESIVE(eng)	19

IŞIK İLE POLİMERİZE OLAN RADYOPAK RESTORATİF NANOHİBRİT ÜNİVERSAL KOMPOZİT DİŞ DOLGUSU

1. ÜRÜN BİLGİSİ

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1. PRODUCT INFORMATION

Universal nano hybrid composite is a filling material that offers high aesthetics, durability and ease of use. Thanks to the combination of micro and nano-sized particles in a hybrid structure, both mechanical properties and aesthetic appearance are optimized. In dentistry, it is suitable for aesthetic restorations on the anterior teeth and functional and durable restorations on the posterior teeth.

2. CONTENT

- Resin Matrix: Bis-GMA, UDMA, TEGDMA based composite resins
- Inorganic Fillers: Nano silica- (0.3-3 µm), quartz.
- Filler Ratio: 80% by weight, 64% by volume
- Photo launcher: Camphorquinone (CQ)
- Pigments
-

3. FEATURES & BENEFITS

- High Aesthetics: Light transmittance compatible with natural tooth structure.
- Mechanical Strength: Resistance to chewing forces.
- Easy Formability: Non-adhesion to tools, fluid and ideal viscosity for modeling.
- Color Options: Natural look with various color tones.
- Long-Term Stability: Abrasion resistance.
- Low Polymerization Shrinkage: Less marginal leakage and long-term stability of tooth-filling fit.

4. INDICATIONS

- In the direct restoration of cavities I-II-III-IV and V.
- In cavities on the surface of the root.
- In the restoration of broken teeth.
- In cases of diastema.
- In direct veneer laminate restoration.
- Core construction.
- Aesthetically enhancing color and shape corrections.

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6. APPLICATION INSTRUCTIONS

- Determine the desired color in daylight with the VITA scale.
- In order for the color choice to be correct, the teeth must be clean.
- Prepare the cavity for composite restorations.
- Isolate the tooth from the teeth next to it. Use transparent thin strips.
- The cavity must be clean and dry. Dry dentin surfaces to strengthen bonding.
- Apply acid and bonding according to the rules.
- Open the lid of the composite filling.
- Apply the composite filling to the cavity in a way that is not thicker than 2 mm.
- Do not remove the oxygen inhibition layer between the incremental layers of the composite to optimize the binding strength to each subsequent layer.
- Keep the beam device as close to the cavity as possible during curing.
- LED light source emitting blue light in a narrow wavelength range between 450-490 nm (1000-1200mW/cm2 power), 20 sec. irradiation for each layer, or HALOGEN light source between 380-500 nm (700mW/cm2 power) 40 sec. irradiation for each layer is recommended as a photo initiator.
- Model the composite in accordance with the natural tooth morphology.
- Put the finishing touches on abrasive discs or burs.
- Complete the polishing process to get a shiny surface.
- Take care to ensure that the airway is free of fat and other contaminants.
- After the application is finished, close the lid of the tube to make sure that it does not get air.
- Closing the lid of the product after each use does not change the current expiration date. The product can be used for the specified shelf life.

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8. PACKAGING

- It is packaged in 2g, 4g, 500g and 1000g.
- It is available in 9 different colors as A1 – A2 – A3 – A3,5 – B1 – B2 -UD – XBW - BW.

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Kids Nanohybrid

KIDS NANO HYBRID COMPOSITE FILLING

LIGHT CURED COLORED RESTORATIVE KIDS NANOHYBRID COMPOSITE DENTAL FILLING MATERIAL

1. PRODUCT INFORMATION

The purpose of the production of colored composites is to facilitate the treatment process in dentistry and to make the experiences of patients (especially children) more positive. Colored nano hybrid composite is a filling material that offers high aesthetics, durability and ease of use. Thanks to the combination of micro and nano-sized particles in a hybrid structure, both mechanical properties and aesthetic appearance are optimized. In dentistry, it is ideal for achieving colorful, aesthetic, functional and durable restorations in milk teeth. Optionally, it is also suitable for use in adult patients.

2. CONTENT

- Resin Matrix: Bis-GMA, UDMA, TEGDMA based composite resins
- Inorganic Fillers: Nano silica- (0.3-3 µm), quartz.
- Filler Ratio: 80% by weight, 64% by volume
- Photo launcher: Camphorquinone (CQ)
- Pigments

3. FEATURES & BENEFITS

- **Reducing Fear of Treatment:** Children are often afraid of dental treatments. Colorful and fun materials are designed to reduce children's fear and make them perceive the treatment as a game.
- **Increasing Motivation:** Since it gives children the opportunity to choose color, they may be more cooperative with the physician during treatment.
- **Visual Difference:** The fact that the colored filling is clearly visible on the tooth facilitates the post-treatment control and follow-up processes. It is easier to remove or replace if necessary.
- **Ideal Mechanical Strength:** Resistance to chewing forces.
- **Easy Formability:** Non-adhesion to tools, fluid and ideal viscosity for modeling.
- **Long-Term Stability:** Abrasion resistance.
- **Low Polymerization Shrinkage:** Less marginal leakage and long-term stability of tooth-filling fit

4. INDICATIONS

- In the restoration of milk teeth.
- In temporary restorations of permanent teeth.
- Fissure sealant applications.
- Optionally, permanent restorations in permanent teeth

5. CONTRAINDICATIONS

- In direct pulp capping cases, in areas very close to the pulp.
- In patients with sensitivity to any of the components of the product.
- Inadequate insulation and humidity control.
- In areas where proper curing cannot be achieved.
- In patients with contact allergies to substances with methacrylate monomer or methacrylate polymer content.

6. APPLICATION INSTRUCTIONS

- Prepare the cavity for composite restorations.
- Isolate the tooth from the teeth next to it. Use transparent thin strips.
- The cavity must be clean and dry. Dry dentin surfaces to strengthen bonding.
- Apply acid and bonding according to the rules.
- Open the lid of the composite filling.
- Apply the composite filling to the cavity in a way that is not thicker than 2 mm.
- Do not remove the oxygen inhibition layer between the incremental layers of the composite to optimize the binding strength to each subsequent layer.
- Keep the beam device as close to the cavity as possible during curing.
- LED light source emitting blue light in a narrow wavelength range between 450-490 nm (1000-1200mW/cm² power), 20 sec. irradiation for each layer, or HALOGEN light source between 380-500 nm (700mW/cm² power) 40 sec. irradiation for each layer is recommended as a photo initiator.
- Model the composite in accordance with the natural tooth morphology.
- Put the finishing touches on abrasive discs or burs.
- Complete the polishing process to get a shiny surface.
- Take care to ensure that the airway is free of fat and other contaminants.
- After the application is finished, close the lid of the tube to make sure that it does not get air.
- Closing the lid of the product after each use does not change the current expiration date. The product can be used for the specified shelf life.

7. ALERTS

- If any nausea or other allergic reactions occur during practice, discontinue use immediately and seek medical attention.
- Uncured methacrylate resin can cause contact dermatitis and damage to the dental pulp. Avoid contact with skin, eyes and soft tissue. If it comes into contact, clean it thoroughly with water.
- Wear personal protective equipment such as disposable masks, gloves and goggles when working with this product.
- If swallowed, drink plenty of water, do not induce vomiting. Consult a doctor immediately.
- Do not use this product for any purpose other than the treatments outlined in the indications section in this user manual.
- This product is intended for use in dentistry only.
- Avoid the use of eugenol, eugenol may disrupt the polymerization of the composite or cause discoloration, so zinc oxide should not be used with cement containing eugenol.
- It may cause discoloration when it comes into contact with cationic mouthwashes or chlorhexidine.
- Clean the application area before use, if it comes into contact with an area you do not want, remove it from the area before polymerization.
- Protect from visible light before use.

8. PACKAGING

- It is packaged in 2g, 4g, 500g and 1000g.
- It is available in Blue, Pink, Purple, Yellow, Green and Orange colors.

9. STORAGE CONDITIONS

- The product should be stored in a cool, dry and dark environment.
- Store syringes between 15 and 25 degrees Celsius at a relative humidity of 30% +/-5%.
- Store dry, cool and not exposed to direct sunlight.
- The expiry date is indicated on all packaging. Do not use after the expiration date.
- Shelf life is 3 years.

10. DISPOSAL INSTRUCTION

Dispose of the products according to local regulations.

11. SYMBOLS

	Production date		Non-sterile mark
	Manufacturer		Store in dry place
	Refer to the introduction for use		It is disposable
	Batch No		Medical device
	Temperature indicator		Dehp free
	Caution		Latex free
	CE Standard		Expiration date
	Protect from sunlight		Relative humidity display

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1. PRODUCT INFORMATION

Nano hybrid posterior composite is a filling material that offers high aesthetics, durability and ease of use. Thanks to the combination of micro and nano-sized particles in a hybrid structure, both mechanical properties and aesthetic appearance are optimized. In dentistry, it is suitable for functional aesthetic and durable restorations on the posterior teeth.

2. CONTENT

- **Resin Matrix:** Bis-GMA, UDMA, TEGDMA based composite resins
- **Inorganic Fillers:** Nano silica- (0.3-3 µm), quartz.
- **Filler Ratio:** 80% by weight, 64% by volume
- **Photo launcher:** Camphorquinone (CQ)
- Pigments
-

3. FEATURES & BENEFITS

- **High Aesthetics:** Light transmittance compatible with natural tooth structure.
- **Mechanical Strength:** Resistance to chewing forces.
- **Easy Formability:** Non-adhesion to tools, fluid and ideal viscosity for modeling.
- **Color Options:** Natural look with various color tones.
- **Long-Term Stability:** Abrasion resistance.
- **Low Polymerization Shrinkage:** Less marginal leakage and long-term stability of tooth-filling fit.

4. INDICATIONS

- In the direct restoration of cavities I-II.
- Core construction.
- In the restoration of worn teeth.
- In the repair of the existing composite filling.
- In erosion and abrasion lesions.

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6. APPLICATION INSTRUCTIONS

- Determine the desired color in daylight with the VITA scale.
- In order for the color choice to be correct, the teeth must be clean.
- Prepare the cavity for composite restorations.
- Isolate the tooth from the teeth next to it. Use transparent thin strips.
- The cavity must be clean and dry. Dry dentin surfaces to strengthen bonding.
- Apply acid and bonding according to the rules.
- Open the lid of the composite filling.
- Apply the composite filling to the cavity in a way that is not thicker than 2 mm.
- Do not remove the oxygen inhibition layer between the incremental layers of the composite to optimize the binding strength to each subsequent layer.
- Keep the beam device as close to the cavity as possible during curing.
- LED light source emitting blue light in a narrow wavelength range between 450-490 nm (1000-1200mW/cm² power), 20 sec. irradiation for each layer, or HALOGEN light source between 380-500 nm (700mW/cm² power) 40 sec. irradiation for each layer is recommended as a photo initiator.
- Model the composite in accordance with the natural tooth morphology.
- Put the finishing touches on abrasive discs or burs.
- Complete the polishing process to get a shiny surface.
- Take care to ensure that the airway is free of fat and other contaminants.
- After the application is finished, close the lid of the tube to make sure that it does not get air.
- Closing the lid of the product after each use does not change the current expiration date. The product can be used for the specified shelf life.

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8. PACKAGING

- It is packaged in 4g, 500g and 1000g.
- It is available in 6 different colors as A1 – A2 – A3 – A3,5 – B1 – B2

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Adres: İvedik OSB Mahallesi 1333. Cadde No:11/2

Yenimahalle/ANKARA - TÜRKİYE TEL: 0(312) 336 36 30

LIGHT CURED FLOWABLE COMPOSITE RESIN

LIGHT CURED RADIOPAQUE FLOW RESTORATIVE COMPOSITE DENTAL FILLING

1. PRODUCT INFORMATION

Flow composite filling material is a low-viscosity, light-curing, radiopaque restorative material. It is specially designed for small cavity restorations, liner (priming), repair procedures and aesthetic corrections in dentistry.

2. CONTENT

The product composition includes methacrylate resins (UDMA, TEGDMA), inorganic filling materials, glass, silica, photo initiator. Inorganic fillers (65% by weight, 52% by volume)

3. FEATURES & BENEFITS

- Flow: With its thin and flow structure, it provides ideal spread in cavities.
- Ideal Thixotropy: Allows controlled filler manipulation
- Radiopacity: It can be easily detected in x-rays, thus providing ease of diagnosis.
- Aesthetic Compatibility: Wide color range that mimics natural tooth color.
- Long-Lasting Durability: Strong bonding feature against micro leakage.
- Easy Application: Provides ease of use with precision tip injectors.

4. INDICATIONS

- In Class III and Class V restorations.
- As a base material in class I and II cavities.
- In the repair of porcelain or composite restorations.
- As a liner under amalgam and composite fillings.
- Together with tooth surface corrections and pit & fissure sealants.

5. CONTRAINDICATIONS

- Restorations bearing large occlusal load.
- Situations that require application in thick layers.
- Routine under gingival applications. In subgingival restorations, the fluid structure of the material can complicate control and affect marginal adaptation.
- Inadequate insulation and humidity control.
- In patients with sensitivity and allergy to any of the direct pulp capping and product ingredients.

6. APPLICATION INSTRUCTIONS

- Determine the desired color in daylight with the VITA scale.
- In order for the color choice to be correct, the teeth must be clean.
- Prepare the cavity for composite restorations.
- Isolate the tooth from the teeth next to it. Use transparent thin strips.
- The cavity must be clean and dry. Dry dentin surfaces to strengthen bonding.
- Apply acid and bonding according to the rules.
- Remove the cap of the injector and insert the tip.
- Apply the composite filling to the cavity in a way that is not thicker than 2 mm.
- Do not remove the oxygen inhibition layer between the incremental layers of the composite to optimize the binding strength to each subsequent layer.
- Keep the light cure device as close to the cavity as possible during curing.
- LED light source emitting blue light in a narrow wavelength range between 450-490 nm (1000-1200mW/cm² power), 20 sec. irradiation for each layer, or HALOGEN light source between 380-500 nm (700mW/cm² power) 40 sec. irradiation for each layer is recommended as a photo initiator.
- Take care to ensure that the airway is free of fat and other contaminants.
- After the application is finished, close the lid of the injector to make sure that it will not get air.
- Closing the lid of the product after each use does not change the current expiration date. The product can be used for the specified shelf life.

7. ALERTS

- If any nausea or other allergic reactions occur during administration, discontinue use immediately and seek medical attention.
- Uncured methacrylate resin can cause contact dermatitis and damage to the dental pulp. Avoid contact with skin, eyes and soft tissue. If it comes into contact, clean it thoroughly with water
- If swallowed, drink plenty of water, do not induce vomiting. Consult a doctor immediately.
- Wear personal protective equipment such as disposable masks, gloves and goggles when working with this product.
- Do not use this product for any purpose other than the treatments outlined in the indications section in this user manual.

- Avoid the use of eugenol, eugenol may disrupt the polymerization of the composite or cause discoloration, so zinc oxide should not be used with cement containing eugenol.
- This product is intended for use in dentistry only.
- To avoid cross-contamination between patients, the application tips are disposable for each patient.
- Do not use damaged syringe tubes.
- Retract the plunger after application to prevent unwanted flow of the composite from the tip.
- Protect from visible light before use.

8. PACKAGING

- Flow composite is presented in the form of 2-3-500-1000 g.
- It is available in 4 different colors as A1-A2-A3-A3,5-.

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Yenimahalle/ANKARA - TÜRKİYE TEL: 0(312) 336 36 30

Bulk Flow

LIGHT CURED BULK FLOW COMPOSITE RESIN

LIGHT CURED RADIOPAQUE BULK FLOW COMPOSITE DENTAL FILLING

1. PRODUCT INFORMATION

Bulk Flow composite is a composite filling material with high flow properties used for restorative procedures in dentistry. It is designed for faster and more effective filling of deep cavities, especially in the posterior regions. Thanks to Bulk fill technology, it allows the material to be used up to 4 mm in a single application layer, which eliminates the need to apply thin layers in conventional composites.

2. CONTENT

The product composition includes methacrylate resins (UDMA, TEGDMA), inorganic filling materials, glass, silica, photo initiator. Inorganic fillers (72.4% by weight, 57% by volume)

3. FEATURES & BENEFITS

- **Flow Structure** : The material easily penetrates the micro-details of the cavities, has high adaptability and minimizes air gaps.
- **High Polymerization Depth**: Thanks to its improved light transmittance, it hardens effectively even in thick layers.
- **Low Shrinkage Stress**: Reduces the stresses that occur during hardening, thus ensuring the longevity of the restoration.
- **Aesthetic Harmony**: It has different color options that adapt to the natural tooth structure.
- **Radiopacity**: It can be clearly seen on X-ray images, which facilitates post-treatment evaluation and follow-up processes.
- **Time Saving**: It reduces the number of layers as it can be polymerized up to 4 mm in a single layer.
- **Ease of Application**: Its flow structure provides comfortable use, especially in narrow or hard-to-reach cavities.
- **Reliable Performance**: Offers low polymerization shrinkage and high mechanical strength.

4. INDICATIONS

- As a base material in class I and II cavities.
- In class V restorations that are not in the aesthetic field.
- As a liner under amalgam and composite fillings.
- Together with tooth surface corrections and pit & fissure sealants.
- In the restoration of the teeth to be crowned.

5. CONTRAINDICATIONS

- Restorations bearing large occlusal load.
- Routine under gingival applications. In subgingival restorations, the fluid structure of the material can complicate control and affect marginal adaptation.
- Inadequate insulation and humidity control.
- In patients with sensitivity and allergy to any of the direct pulp capping and product ingredients.

6. APPLICATION INSTRUCTIONS

- Determine the desired color in daylight with the VITA scale.
- In order for the color choice to be correct, the teeth must be clean.
- Prepare the cavity for composite restorations.
- Isolate the tooth from the teeth next to it. Use transparent thin strips.
- The cavity must be clean and dry. Dry dentin surfaces to strengthen bonding.
- Apply acid and bonding according to the rules.
- Open the lid of the composite filling.
- Apply the composite filling to the cavity in a way that is not thicker than 4 mm.
- Do not remove the oxygen inhibition layer between the incremental layers of the composite to optimize the binding strength to each subsequent layer.
- Keep the light cure device as close to the cavity as possible during curing.
- LED light source emitting blue light in a narrow wavelength range between 450-490 nm (1000-1200mW/cm² power), 20 sec. irradiation for each layer, or HALOGEN light source between 380-500 nm (700mW/cm² power) 40 sec. irradiation for each layer is recommended as a photo initiator.
- Take care to ensure that the airway is free of fat and other contaminants.
- After the application is finished, close the lid of the injector to make sure that it will not get air.
- Closing the lid of the product after each use does not change the current expiration date. The product can be used for the specified shelf life.
- For repeated use, make sure that the lid is closed correctly after use.

7. ALERTS

- If any nausea or other allergic reactions occur during administration, discontinue use immediately and seek medical attention.
- Uncured methacrylate resin can cause contact dermatitis and damage to the dental pulp. Avoid contact with skin, eyes and soft tissue. If it comes into contact, clean it thoroughly with water
- If swallowed, drink plenty of water, do not induce vomiting. Consult a doctor immediately.
- Wear personal protective equipment such as disposable masks, gloves and goggles when working with this product.

- Do not use this product for any purpose other than the treatments outlined in the indications section in this user manual.
- Avoid the use of eugenol, eugenol may disrupt the polymerization of the composite or cause discoloration, so zinc oxide should not be used with cement containing eugenol.
- This product is intended for use in dentistry only.
- To avoid cross-contamination between patients, the application tips are disposable for each patient.
- Do not use damaged syringe tubes.
- Retract the plunger after application to prevent unwanted flow of the composite from the tip.
- Protect from visible light before use.

8. PACKAGING

- Bulk Flow composite is offered in the 2g. syringe tubes.
- It is available in 3 different colors as A1–A2–CLEAR.

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Flor Liner

LIGHT CURED LINING MATERIAL

RADIOPAQUE FLOWABLE RESIN MODIFIED GLASS IONOMER FILLING AND LINING MATERIAL

1. PRODUCT INFORMATION

Flor releasing resin modified glass is an ionomer-based filling and base material. It is a low-viscosity, light-curing, radiopaque restorative material.

2. CONTENT

In the product composition; methacrylate resins (UDMA, TEGDMA), inorganic filling materials, fluorosilicate glass, quartz, silica, photo initiator and inhibitor. Inorganic fillers (65% by weight, 52% by volume)

3. FEATURES & BENEFITS

- **Flor Release:** Protects the tooth against caries.
- **Flow:** With its thin and flow structure, it provides ideal spread in cavities.
- **Ideal Thixotropy:** Allows controlled filler manipulation
- **Rodiopacity:** It can be easily detected in x-rays, thus providing ease of diagnosis.
- **Easy Application:** Provides ease of use with precision tip injectors

4. INDICATIONS

- Base material under amalgam and composites.
- Pit & Fissure sealant.
- Liner in deep cavities.
- Treatment of small lesions.

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6. APPLICATION INSTRUCTIONS

- Prepare the cavity for restorations.
- Isolate the tooth from the teeth next to it. Use transparent thin strips.
- A calcium hydroxide-based base material should be used in areas close to the pulp.
- The cavity must be clean and dry. Dry dentin surfaces to strengthen bonding.
- Apply acid and bonding according to the rules.
- Remove the cap of the injector and insert the tip.
- Apply the restorative material to the cavity no thicker than 1 mm.
- Keep the beam device as close to the cavity as possible during curing.
- LED light source emitting blue light in a narrow wavelength range between 450-490 nm (1000-1200mW/cm² power), 20 sec. irradiation for each layer, or HALOGEN light source between 380-500 nm (700mW/cm² power) 40 sec. irradiation for each layer is recommended as a photo initiator.
- Take care to ensure that the airway is free of fat and other contaminants.
- After the application is finished, close the lid of the injector to make sure that it will not get air.
- Closing the lid of the product after each use does not change the current expiration date. The product can be used for the specified shelf life.

7. ALERTS

- If any nausea or other allergic reactions occur during application, discontinue use immediately and seek medical attention.
- Uncured methacrylate resin can cause contact dermatitis and damage to the pulp. Avoid contact with skin, eyes and soft tissue. If it comes into contact, clean it thoroughly with water
- If swallowed, drink plenty of water, do not induce vomiting. Consult a doctor immediately.
- Wear personal protective equipment such as disposable masks, gloves and goggles when working with this product.
- Do not use this product for any purpose other than the treatments outlined in the indications section in this user manual.
- Avoid the use of eugenol, eugenol may disrupt the polymerization of the restorative or cause discoloration, so zinc oxide should not be used with cements containing eugenol.
- This product is intended for use in dentistry only.
- To avoid cross-contamination between patients, the application tips are disposable for each patient.
- Do not use damaged syringe tubes.

- Retract the plunger after application to prevent unwanted flow of the restorative from the tip.
- Protect from visible light before use.

8. PACKAGING

Flor liner: Available in 2 g or 4 x 2 g syringe tubes.

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Yenimahalle/ANKARA - TÜRKİYE TEL: 0(312) 336 36 30

Orthodontic Flow

LIGHT CURED NANO HYBRID ORTHODONTIC ADHESIVE

LIGHT CURED COMPOSITE-BASED ORTHODONTIC BRACKET ADHESIVE

1. PRODUCT INFORMATION

It is an orthodontic dental material developed on the basis of composite, polymerized by light, with high adhesion strength and durability, used to fix orthodontic brackets on the tooth surface. Its viscosity is adjusted in such a way that it can easily spread to the pores of the adhesion surface of the bracket and prevent slipping during tooth placement. It is compatible with both metal and ceramic brackets and can be used safely in long-term orthodontic treatments.

2. CONTENT

The product composition includes methacrylate resins (UDMA, TEGDMA), inorganic filling materials, glass, silica, photoinitiator. (70% by weight, 56% by volume)

3. FEATURES & BENEFITS

- **High Adhesion Strength:** Creates a strong bond between the tooth and the bracket.
- **Controlled Working Time:** Since it is polymerized with light, it gives the physician sufficient working time.
- **Ideal Viscosity:** The adhesion surface of the bracket is adjusted so that it can easily spread into the pores and prevent slipping during tooth placement.
- **Low Polymerization Shrinkage:** Thanks to the high filler ratio it contains, it shows minimum shrinkage effect during hardening.
- Compatible with routine bracket bonding protocols.
- It is compatible with natural tooth color.
- It can be removed without damaging the enamel of the tooth.
- It is one-component, does not require mixing.

4. INDICATIONS

- **Fixing of Orthodontic Brackets:** Bonding metal, ceramic, composite and polycarbonate brackets to the tooth surface.
- **Bonding Orthodontic Tapes:** Fixing special accessories or tubes.
- Fixation of temporary orthodontic devices.

5. CONTRAINDICATIONS

- Methacrylate is contraindicated in patients with a history of severe allergic reactions to resins.
- **Unsuitable Tooth Surfaces:** Weak enamel, dentin surface, and decayed or plaque-covered tooth surfaces.
- **Inadequate Insulation:** Situations where humidity control cannot be achieved.

6. APPLICATION INSTRUCTIONS

- **Preparation of the Tooth Surface:**
- **Cleaning:** Tooth surfaces are cleaned with a polishing brush and pumice. This removes plaque and organic waste from the surface. After cleaning, the surface is dried with water and air spray.
- **Insulation:** Insulation of the tooth surface is provided with the help of rubber cover or cotton rolls. Insulation increases adhesion success by providing moisture control.
- **Etching:** 35-37% phosphoric acid is applied to the enamel surface of the tooth. The acid is left on the enamel surface **for 15-30 seconds**. After the application of acid, the surface is washed with plenty of water and slightly dried. The surface should take on a matte appearance.
- Bonding process:**
- **Bonding Agent Application:** The bonding agent is applied as a thin layer to the acidified tooth surface. With a mild air spray, the bonding agent is diluted and the solvent is evaporated. An LED light source (1000-1200mW/cm2 power) emitting blue light in a narrow wavelength range between 450-490 nm, 10 sec. irradiation for each tooth, or a HALOGEN light source between 380-500 nm (700mW/cm2 power) for 20 seconds for each layer is recommended
- **Application of Adhesive to Bracket:** Composite-based adhesive is applied in an even layer to the bracket base (mesh surface). Care is taken not to overdo the adhesive.
- **Placement of the Bracket:** The adhesive bracket is placed on the tooth surface in the appropriate position. By applying light pressure, the excess adhesive is carried away from the edges.
- **Cleaning the Overflowing Adhesive:** Excess is removed with the help of a fine probe or micro brush.
- **Curing:** A LED light source (1000-1200mW/cm2 power) emitting blue light in a narrow wavelength range between 450-490 nm, 10 sec. irradiation for each tooth, or a HALOGEN light source between 380-500 nm (700mW/cm2 power) for 20 seconds for each layer is recommended

7. ALERTS

- If any nausea or other allergic reactions occur during application, discontinue use immediately and seek medical attention.
- Uncured methacrylate resin can cause contact dermatitis and damage to the pulp. Avoid contact with skin, eyes and soft tissue. If it comes into contact, clean it thoroughly with water
- If swallowed, drink plenty of water, do not induce vomiting. Consult a doctor immediately.
- Wear personal protective equipment such as disposable masks, gloves and goggles when working with this product.
- Do not use this product for any purpose other than the treatments outlined in the indications section in this user manual.
- Avoid the use of eugenol, eugenol may disrupt the polymerization of the restorative or cause discoloration, so zinc oxide should not be used with cements containing eugenol.
- This product is intended for use in dentistry only.
- Attention should be paid to the rules of disinfection so as not to lead to cross-contamination between patients
- Do not use damaged syringe tubes.
- Retract the plunger after application to prevent unwanted flow of the restorative from the tip.
- Protect from visible light before use.

8. PACKAGING

Orthodontic Flow composite is packaged as 2g, 3g, 500g and 1000g.

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Relative humidity display

Warning
H3 15
H335
H3 19

It can lead to skin, eye and respiratory tract irritation.

12. COMPANY INFORMATION

Manufacturer: FSM DENTAL MEDIKAL MAKINA SANAYI VE TIC.LTD.

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Yenimahalle/ANKARA - TÜRKİYE TEL: 0(312) 336 36 30

Bond

LIGHT CURED 7th. GENERATION DENTAL ADHESIVE

LIGHT CURED ONE-COMPONENT 7th. GENERATION UNIVERSAL DENTAL ADHESIVE

1. PRODUCT INFORMATION

The 7th generation one-component bond system is a self-etch bonding system in dentistry in which both pickling and primary and bonding agents are combined. These systems provide effective bonding to both dentin and enamel and offer ease of application and time savings.

It is suitable for use in selective etch (roughening only tooth enamel with phosphoric acid) and total etch (roughening on both enamel and dentin with phosphoric acid) procedures. It is a universal bonding agent that polymerizes with radiation with ideal adhesion and adhesion strength.

2. CONTENT

UDMA, TEDGMA, 4-META, 10-MDP, Camphorquinone, acetone, silane-silica-water.

3. FEATURES & BENEFITS

- **Easy and Fast Application:** Since it is a single component, it shortens the application time.
- **Minimal Technical Precision:** Humidity control requires less precision than conventional total-etch systems.
- **Reduced Post-Operative Sensitivity:** It provides obstruction of dentin tubules and reduces the risk of sensitivity.
- **Low Film Thickness:** The thin film layer increases the adaptability of the restoration.
- **High Bonding Strength:** Provides strong bonding to both enamel and dentin.

4. INDICATIONS

- Bonding of direct restorations of all cavity classes (Black) using light, dual or self-polymerizing methacrylate-based composites/compomers.
- Bonding of light, dual or self-polymerizing core making materials.
- Covering hypersensitive tooth areas.
- Covering cavities before amalgam restorations.
- Bonding fissure sealants.
- Closure of cavities and core preparations before applying temporary cementation to indirect restorations (in accordance with the dentin closure technique applied without delay)
- Cementation of indirect restorations with light, dual or self-polymerizing adhesive resin cements.
- Intraoral repair of composite and composite restorations, metal-supported porcelain and all-ceramic restorations, metal restorations.

5. CONTRAINDICATIONS

- Use of this product in cases of known or suspected allergy to the ingredients of this product.
- In case of known or suspected allergy to methacrylate components.
- Do not allow Bond to come into direct contact with unprotected pulp tissue

6. APPLICATION INSTRUCTIONS

CAVITY PREPARATION

- Complete tooth decay cleaning and cavity shaping.
- Clean the cavity with an air-water spray and leave it slightly damp (do not dry it completely). **Over drying can lead to dehydration of the dentin tubules and a decrease in their bonding strength.**

ISOLATION

- Rubber dam: Recommended for humidity control.
- Alternatively, cotton rolls and saliva absorbers can also be used.

BOND APPLICATION

- **Preparation of the Bond agent:** Shake the bottle well to make it homogeneous.
- **Application of Bond Agent:** Apply the bond directly to the enamel and dentin surface with a micro brush or applicator.
- Gently feed the Bond agent into the surface and apply it with scrubbing motions for about 20 seconds.
- **Light Air Ventilation:** Gently dispense the Bond with low-pressure air, forming a thin layer on the surface. Do not use excessive air pressure; This can cause Bond to move away from the surface.
- **Curing:** An LED light source (1000-1200mW/cm² power) emitting blue light in a narrow wavelength range between 450-490 nm is recommended as a photo initiator, 10 sec. irradiation for each layer, or a HALOGEN light source between 380-500 nm (700mW/cm² power) for 20 sec. irradiation for each layer.

COMPOSITE APPLICATION

- Apply your composite restoration directly on the cured bond layer.
- Polymerize the composite by layering in accordance with the manufacturer's instructions

7. ALERTS

- Since the bond polymerizes with light, it must be protected against ambient light. After use, the lid should be closed immediately.
- Materials containing eugenol can adversely affect the polymerization of bond.
- If a problem is encountered in allergic and hypersensitive people, the application should be stopped immediately, washed with plenty of water and a specialist doctor should be consulted.
- Personal protective equipment such as gloves, face masks, and safety goggles should be worn at all times.
- Avoid inhalation or ingestion of material.
- Bond is flammable and volatile. Do not use directly next to flammable. Keep away from explosive sources. Keep away from direct sunlight.
- In case of contact with eyes, rinse immediately with plenty of water or seek medical attention.
- It is not recommended to use disinfectants, hemostats and rinse solutions containing EDTA or hydrogen peroxide. Residues of these substances can adversely affect adhesion and adhesive polymerization reaction. Chlorhexidine solutions can be used for disinfection.
- Dispose of the contents and bottle according to local regulations.

8. PACKAGING

.It is packaged in 2ml, 5ml, 500ml and 1000ml.

9. STORAGE CONDITIONS

- Store syringes between 15° and 25 ° Cat +/-5 30% rh.
- Do not use damaged syringe tubes.
- If it will not be used for a long time, store it in a cool and dark environment.
- The shelf life of the product is 2 years

10. SYMBOLS

Production date

Non-sterile mark

Manufacturer

Store in dry place

Refer to the introduction for use

It is disposable

Batch No

Medical device

Temperature indicator

Dehp free

Caution

Latex free

CE Standard

Expiration date

Protect from sunlight

Relative humidity display

11. COMPANY INFORMATION

Manufacturer: FSM DENTAL MEDİKAL MAKİNA SANAYİ VE TİC.LTD.
Address: İvedik OSB Mahallesi 1333. Cadde No:11/2
Yenimahalle/ANKARA,TURKIYE
Phone: 0(312) 336 36 30
Web Address: www.fsmdental.com